

Balance Control

Topics Covered


- Control enabling logic.
- PID-based balance control.

Prerequisites


- Filtering laboratory experiment.
- PD Control laboratory experiment.
- Pendulum Moment of Inertia laboratory experiment.
- Rotary pendulum module is attached to the QUBE-Servo 2.

2 In-Lab Exercises

Based on the VI developed in the Pendulum Moment of Inertia laboratory experiment, build a VI similar to that shown in Figure 2.1 that balances the pendulum on the QUBE-Servo 2 rotary pendulum using the PD control detailed in the background section of this lab.



(a) Front Panel



(b) Block Diagram

Figure 2.1: VI to run PD balance controller

- Using the VI you made in the Pendulum Moment of Inertia laboratory experiment, construct the controller shown in Figure 2.1:
 - The Counts to Angles subsystem contains the same blocks used in the Pendulum Moment of Inertia laboratory experiment to convert encoder counts to radians. Make sure you use the inverted pendulum angle.
 - To find the velocity of the rotary arm and pendulum, add the high-pass filters $100s/(s + 100)$ similarly to the Filtering laboratory experiment.

- Add the necessary Sum and Gain blocks to implement the PD control given in Equation 1.1.
 - The controller should only be enabled when the pendulum is $\pm 10^\circ$ about the upright vertical position (or ± 0.175 rad). Add Absolute Value, Constant Comparison and Selector blocks to implement this.
2. Set the PD gains as follows: $k_{p,\theta} = -2$ V/rad, $k_{p,\alpha} = 30$ V/rad, $k_{d,\theta} = -2$ V/(rad/s), and $k_{d,\alpha} = 2.5$ V/(rad/s).
 3. Make sure the pendulum is not perturbed in the downward position before starting the controller. Run the VI.
 4. Manually rotate the pendulum in the upright position until the controller engages. The waveform charts should read something similar as shown in Figure 2.2. Attach response of the rotary arm, pendulum, and controller voltage.




Figure 2.2: QUBE-Servo 2 rotary pendulum response

5. As the pendulum is being balanced, describe the responses in the *Rotary Arm (deg)* and the *Pendulum Angle (deg)* waveform chart.
6. Locate the Summation block that precedes `kp_theta`. Vary the Constant block that is connected to the positive input of the Sum block. **Do not set the value too high, keep it within ± 45 to start.** Observe the response in the *Arm Angle (deg)* waveform chart. What variable does this represent in the balance control?
7. Click on the Stop button to stop the VI.

© 2017 Quanser Inc., All rights reserved.

Quanser Inc.
119 Spy Court
Markham, Ontario
L3R 5H6
Canada
info@quanser.com
Phone: 1-905-940-3575
Fax: 1-905-940-3576

Printed in Markham, Ontario.

For more information on the solutions Quanser Inc. offers, please visit the web site at:
<http://www.quanser.com>

This document and the software described in it are provided subject to a license agreement. Neither the software nor this document may be used or copied except as specified under the terms of that license agreement. Quanser Inc. grants the following rights: a) The right to reproduce the work, to incorporate the work into one or more collections, and to reproduce the work as incorporated in the collections, b) to create and reproduce adaptations provided reasonable steps are taken to clearly identify the changes that were made to the original work, c) to distribute and publically perform the work including as incorporated in collections, and d) to distribute and publicly perform adaptations. The above rights may be exercised in all media and formats whether now known or hereafter devised. These rights are granted subject to and limited by the following restrictions: a) You may not exercise any of the rights granted to You in above in any manner that is primarily intended for or directed toward commercial advantage or private monetary compensation, and b) You must keep intact all copyright notices for the Work and provide the name Quanser Inc. for attribution. These restrictions may not be waved without express prior written permission of Quanser Inc.